# Electronic Pulse Stimulator (TENS Unit)

**Operating Manual  
PL-029U Edition V1.0**

# Table of Contents

Introduction.....	1
Indications for Use.....	1
Safety Warning.....	2
Contraindications.....	2
Warnings.....	2
Precautions.....	2
Adverse Reactions.....	4
Environmental Condition for Transport and Storage.....	4
Symbols interpretation.....	5
Safety Test Standards.....	5
Electromagnetic Compatibility and FCC Compliance statement.....	6
Product Specifications.....	12
Setup.....	15
Operating Instruction.....	16
Recommended practice.....	17
Cleaning and Maintenance.....	18
Trouble Shooting.....	19
Recommended Use Positions.....	20
Contact Information.....	21
Électrostimulateur (Appareil TENS).....	22
Elektronischer Pulsstimulator (TENS Gerät).....	27
Stimolatore Elettronico a Impulsi (Unità TENS).....	31
Estimulador eléctrico transcutáneo (Unidad de TENS).....	36

## Contents of Different Language Versions

English.....	1-21
French.....	22-36
German.....	27-30
Italian.....	31-35
Spanish.....	36-40

## **Introduction**

---

Electronic Pulse Stimulator (TENS) delivers electric impulses to tired and sore muscles. These different frequencies of impulses covering Transcutaneous Electrical Nerve Stimulation mimic the action potential coming from the central nervous system to trigger contraction of the muscle. It may be helpful in relieving aches and pains in various parts of the body such as the waist, shoulders, joints, hands and feet.

## **Indications for Use**

To be used for temporary relief of pain associated with sore and aching muscles in the shoulder, waist, back, arm, and leg, due to strain from exercise or normal household and work activities.

## **Safety Warning**

---

### **Contraindications**

Do not use this device on patients who have a cardiac pacemaker, implanted defibrillator, or other implanted metallic or electronic device, because this may cause electric shock, burns, electrical interference, or death.

Do not use this device on patients whose pain syndromes are undiagnosed.

### **Warnings**

Do not apply stimulation over the patient's neck because this could cause severe muscle spasms resulting in closure of the airway, difficulty in breathing, or adverse effects on heart rhythm or blood pressure.

Do not apply stimulation across the patient's chest, because the introduction of electrical current into the chest may cause rhythm disturbances to the patient's heart, which could be lethal.

Do not apply stimulation over, or in proximity to, cancerous lesions.

Do not apply stimulation when the patient is in the bath or shower.

If you have one of the following conditions, please consult with your physician before purchasing or using this device.

Acute disease, malignant tumor, infective disease, pregnant, heart disease, high fever, abnormal blood pressure, lack of skin sensation or an abnormal skin condition, any condition requiring the active supervision of a physician.

### **Precautions**

Do not use this device while driving.

Do not use this device while sleeping.

Do not use this device in high humidity areas such as a bathroom.

Keep the device away from wet, high temperature and direct-sunlight place.

Keep this device out of reach of children.

Stop using this device at once if you feel pain, discomfort, dizziness or nausea and consult your physician.

Do not attempt to move the electrode pads while the device is operating.

Do not use the device around the heart, on the head, mouth, pudendum or blemished skin areas.

Do not apply stimulation of this device in the following conditions:

(1) across the chest because the introduction of electrical current into the chest may cause rhythm disturbances to the heart, which could be lethal;

(2) over painful areas. Please consult with your physician before using this device if you have painful areas;

(3) over open wounds or rashes, or over swollen, red, infected, or inflamed areas or skin eruptions (e.g.,

phlebitis, thrombophlebitis, varicose veins). Apply stimulation only to normal, intact, clean, healthy skin;

(4) in the presence of electronic monitoring equipment (e.g., cardiac monitors, ECG alarms). The electronic stimulator may not operate properly when the electrical stimulation device is in use;

(5) while operating machinery, or during any activity in which electrical stimulation can put you at risk of injury;

(6) on children.

Be aware of the following.

(1) to consult with your physician before using this device. The simulation with the device may:

i. cause lethal rhythm disturbances to the heart in susceptible individuals;

ii. disrupt the healing process after a recent surgical procedure;

(2) that the device is not effective for pain of central origin, including headache;

(3) that the device is not a substitute for pain medications and other pain management therapies;

(4) that the device has no curative value;

(5) that the device is a symptomatic treatment and, as such, suppresses the sensation of pain that would otherwise serve as a protective mechanism;

(6) that the long-term effects of electrical stimulation are unknown;

(7) that the user may experience skin irritation, burns or hypersensitivity due to the electrical stimulation or electrical conductive medium (gel);

(8) if the user has suspected or diagnosed epilepsy, the user should follow precautions recommended by his or her physician;

(9) to use caution if the user has a tendency to bleed internally, such as following an injury or fracture;

(10) use caution if stimulation is applied over the menstruating uterus;

(11) use caution if stimulation is applied over areas of skin that lack normal sensation;

(12) stop using the device if the device does not provide pain relief;

(13) use this device only with the leads, electrodes, and accessories that the manufacturer recommends.

(14) Do not share the use of the electrode pads with others.

(15) Do not use the device while it's charging.

(16) The device contains the lithium battery. If overheating of the device occurred during the charging, stop the charging or operation immediately and report to the seller.

(17) Dispose of the battery-containing device according to the local, state, or federal laws.

The long-term effects of electrical stimulation are unknown.

Since the effects of stimulation of the brain are unknown, stimulation should not be applied across the head, and electrodes should not be placed on opposite sides of the head.

The safety of electrical stimulation during pregnancy has not been established.

Some patients may experience skin irritation or hypersensitivity due to the electrical stimulation or

electrical conductive medium (gel).

Patients with suspected or diagnosed heart disease should follow precautions recommended by their physicians.

Patients with suspected or diagnosed epilepsy should follow precautions recommended by their physicians.

Use caution if stimulation is applied over the menstruating or pregnant uterus.

### **Adverse Reactions**

Patients may experience skin irritation and burns beneath the stimulation electrodes applied to the skin;

Patients may experience headache and other painful sensations during or following the application of electrical stimulation near the eyes and to the head and face.

Patients should stop using the device and should consult with their physicians if they experience adverse reactions from the device.

### **Environmental condition for normal working, transport and storage**

- Normal working ambient temperature: 5~40°C
- Normal working ambient humidity: 15%~90% RH
- Store and transport ambient temperature: -25 ~70°C
- Store and transport ambient humidity: 0%~90% RH
- Atmospheric pressure: 70~106kPa

## Symbols interpretation

	Fragile, handle with care		Type BF applied part
	Keep the product in the dry place Away from water and rain.		CAUTION, Avoid injury. Read and understand owner's manual before operating this product.
	This way up		Manufacturer
	Product package should be recycled.		Symbol for "AUTHORISED REPRESENTATIVE IN THE EUROPEAN COMMUNITY"
	Unrecyclable		Batch code
	Date of manufacture	FDA 510(k)	FDA 510(k) cleared
	Serial number	IP22	IP code of the device

### Safety Test Standards:

- Medical Devices Directive 93/42/EEC
- IEC 60601-1:2005+A1:2012/EN 60601-1:2006 Medical electrical equipment - Part 1: General requirements for basic safety and essential performance
- IEC 60601-1-2:2007/EN 60601-1-2:2007 Medical electrical equipment - Part 1-2: General requirements for safety - Collateral standard: Electromagnetic compatibility - Requirements and tests
- IEC 60601-2-10:2012/EN 60601-2-10:2000+A1:2001 Medical electrical equipment - Part 2-10: Particular requirements for the safety of nerve and muscle stimulators
- IEC 60601-1-11:2010 Medical electrical equipment -- Part 1-11: General requirements for basic safety and essential performance -- Collateral standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment.
- EN 980 Symbols for use in the labeling of medical devices
- EN 1041 Information supplied by the manufacturer with medical devices
- IEC/60601-1-6/ EN 60601-1-6 Medical electrical equipment – Part1-6: General requirements for basic safety and essential performance – Collateral standard: Usability

- IEC 60601-1-11/ EN 60601-1-11 Medical electrical equipment – Part 1-11: General requirements for basic safety and essential performance – Collateral standard: Requirements for medical electrical equipment and medical electrical systems used in home healthcare environment
- IEC 62304/ EN 62304 Medical device software - Software life-cycle processes
- IEC 62366/ EN 62366 Medical devices – Application of usability engineering to medical devices
- ISO 10993-1 Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process

### **Electromagnetic Compatibility and FCC Compliance Statement**

- 1) This product needs special precautions regarding electromagnetic compatibility (EMC) and needs to be installed and put into service according to the EMC information provided, and this unit can be affected by portable and mobile radio frequency (RF) communications equipment.
- 2) Do not use a mobile phone or other devices that emit electromagnetic fields, near the unit. This may result in incorrect operation of the unit.
- 3) Caution: This unit has been thoroughly tested and inspected to assure proper performance and operation!
- 4) Caution: This machine should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, this machine should be observed to verify normal operation in the configuration in which it will be used

<b>Guidance and manufacture’s declaration – electromagnetic emission</b>		
The device <i>is</i> intended for use in the electromagnetic environment specified below. The customer of the user of the device should assure that it is used in such an environment.		
<b>Emission test</b>	<b>Compliance</b>	<b>Electromagnetic environment – guidance</b>
RF emissions CISPR 11	Group 1	The device use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emission CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable (internal battery powered)	

Voltage fluctuations/ flicker emissions  IEC 61000-3-3	Not applicable  (internal battery powered)	
---	--	--

<b>Guidance and manufacture's declaration – electromagnetic immunity</b>			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of device should assure that it is used in such an environment.			
<b>Immunity test</b>	<b>IEC 60601 test level</b>	<b>Compliance level</b>	<b>Electromagnetic environment - guidance</b>
Electrostatic discharge (ESD)  IEC 61000-4-2	±6 kV contact  ±8 kV air	±6 kV contact  ±8 kV air	Floors should be wood, concrete or ceramic tile. If floor are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst  IEC 61000-4-4	±2 kV for power supply lines  ±1 kV for input/output lines	Not applicable  (internal battery powered)	Mains power quality should be that of a typical commercial or hospital environment.
Surge  IEC 61000-4-5	± 1 kV line(s) to line(s)  ± 2 kV line(s) to earth	Not applicable  (internal battery powered)	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines  IEC 61000-4-11	<5% $U_T$  (>95% dip in $U_T$ )  for 0.5 cycle  40% $U_T$  (60% dip in $U_T$ )  for 5 cycles  70% $U_T$	Not applicable  (internal battery powered)	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.

	(30% dip in $U_T$ ) for 25 cycles  <5% $U_T$ (>95% dip in $U_T$ ) for 5 sec		
Power frequency (50Hz/60Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE $U_T$ is the a.c. mains voltage prior to application of the test level.			

<b>Guidance and manufacture's declaration – electromagnetic immunity</b>			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
<b>Immunity test</b>	<b>IEC 60601 test level</b>	<b>Compliance level</b>	<b>Electromagnetic environment - guidance</b>
Conducted RF IEC 61000-4-6	3 V <sub>rms</sub> 150 kHz to 80 MHz	3 V <sub>rms</sub>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p><b>Recommended separation distance</b></p> $d = 1,2\sqrt{P}$

<p>Radiated RF IEC 61000-4-3</p>	<p>3 V/m 80 MHz to 2.5 GHz</p>	<p>3 V/m</p>	<p> <math>d = 1,2\sqrt{P}</math> 80 MHz to 800 MHz  <math>d = 2,3\sqrt{P}</math> 800 MHz to 2,5 GHz </p> <p>Where <math>P</math> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <math>d</math> is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,<sup>a</sup> should be less than the compliance level in each frequency range.<sup>b</sup></p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
--------------------------------------	------------------------------------	--------------	---

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.

b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

**Recommended separation distances between  
portable and mobile RF communications equipment and the device .**

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 KHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

The subject device has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.

The product generates, uses, and can radiate radio frequency energy and, if not installed and used accordance with the instructions, may cause harmful interference to radio communications.

However, there is no guarantee that the interference will not occur in a particular installation. If the product does cause harmful interference to radio or television reception, which can be determined by turning the product on or off, the user is encouraged to try to correct the interference by one or more of the following measures:

- a) Reorient or relocate the receiving antenna;
- b) Increase the separation between the product and the receiver;
- c) Consult the dealer or an experienced radio / TV technician for help.
- d) Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.

Changes or modifications to this product not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

## Product Specifications

---

### Accessories included in the package.

- (1).Tens unit controller \* 1pc
- (2).Large gel pad\*2pcs; palm shape gel pad\*4pcs; small gel pads \* 2pcs
- (3).Output wire \* 3pcs
- (4).USB line \*1pc
- (5).Gel pad collector \* 1pc
- (6).Manual \*1pc
- (7).Collecting bag\*1pc

### Technical Information

Model/type	PL-029U	Weight	110g
Power supply	Powered by internal 3.7V li-ion battery	Automatic shutoff	20 minutes
Waveform and wave shape	Biphasic rectangular wave pulse	Degree of protection against electric shock	Type BF applied part
Pulse duration	100us (Microseconds)	Type of protection against electric shock	Internally powered equipment
Pulse frequency	1-100Hz (Hz=vibration per second)	Grade of waterproof	IP22
Output Voltage	Max. 50pp $\pm$ 20%(at 500ohm load)	Product life	1 year
Treatment time	10,20,30,40,50,60minutes	Lifetime for electrode	Storage for 2 year (no use) , Times of reusable: 30 times
Output intensity	0 to 20 levels, adjustable	Mode of operation	Continuous operation
Modes	18 auto modes	Software version	A0
Typical operation time of Battery	If to use both channels at level 20, the battery can be used for about 180 minutes after fully charged.	The time required for me equipment to warm from the minimum storage temperature between uses until it is ready for	30 minutes

		intended use	
Behaviour of me equipment while the rechargeable internal electrical power source is charging:	The battery icon on the right corner will flash during charging and will be still with full capacity after fully charged.	The time required for me equipment to cool from the maximum storage temperature between uses until it is ready for intended use	15 minutes
Typical service life of Battery	300 times of recharging	Adapter for charging	Please use output DC5V and output current 0.3-2.0A adapter for charging
Note: Not intended to be sterilized.			
Not for use in an OXYGEN RICH ENVIRONMENT			


## Programs

Program name	Time min.	Frequency (Hz)	Pulse Width ( $\mu$ s)
Mode 1	10,20,30,40,50,60	10-35	100
Mode 2	10,20,30,40,50,60	8.75-60	100
Mode 3	10,20,30,40,50,60	4-100	100
Mode 4	10,20,30,40,50,60	18-100	100
Mode 5	10,20,30,40,50,60	18-52.5	100
Mode 6	10,20,30,40,50,60	1-100	100
Mode 7	10,20,30,40,50,60	1-35	100
Mode 8	10,20,30,40,50,60	6.6-10	100
Mode 9	10,20,30,40,50,60	4-8	100
Mode 10	10,20,30,40,50,60	62.5	100
Mode 11	10,20,30,40,50,60	2	100
Mode 12	10,20,30,40,50,60	12.5-55.5	100
Mode 13	10,20,30,40,50,60	100	100
Mode 14	10,20,30,40,50,60	80	100
Mode 15	10,20,30,40,50,60	50	100
Mode 16	10,20,30,40,50,60	6	100
Mode 17	10,20,30,40,50,60	30	100
Mode 18	10,20,30,40,50,60	Combination of the above	100

## Setup

Unpack the box of the product, take the product and accessories out, and snap the connect unit onto the electrode pad.

## Operating instruction

---

The following steps are used to guide the device operation.

- Electronic Pulse Stimulator (TENS) needs to be charged for up to 8 hours before the first use.
- Connect a pair of electrode pads to one connecting wire by snapping them on; the other end of the connecting wire is connected to the left output of the device. Similarly, the other pair of electrode pads is connected to the remaining connecting wire and the right output of the device.
- Attach one pair or two pair of the electrode pads to the treatment area, such as shoulder and legs.
- Press the ON/OFF button for 2 seconds to turn the unit on, the unit will automatically start at Mode 1 and Channel A.
- Lightly touching the screen to change mode as your like, make sure not change the mode too fast.
- Before changing intensity of mode, firstly you should know which channel you are at. Since the intensity changing of channel A and channel B is separate. If icon A is flashing, you are can change the intensity of channel A only. If icon B is flashing, you can change the intensity of channel B only.
- If you want to change from channel A from channel B, you can press the button A/B once or you can just simply press A or B icon on the two sides of intensity bars at bottom of the screen to change the channel.
- When certain mode has been chosen, you can change the intensity of this mode. The icon A will be flashing, you can change the intensity of channel A now. If you want to change the intensity of channel B, press the button A/B once or you can just simply press A or B icon, then the icon B is flashing, now you can change the intensity of channel B only.
- There are two ways to change the intensity, you can touch the icon + on the screen OR pressing the + button on the right side of uni to increase the intensity. By touching the icon- on the screen or press the - button on the right side of the unit, you can decrease the intensity.
- When pressing the SUSPEND button, the unit will pause and there is not output till you press the button SUSPEND again. In SUSPEND status, the intensity, mode, timer etc are not adjustable. But the screen lock function is available.
- If you like certain mode and you want to keep using it, you can simply press the cycle icon  on top of the screen, then the current mode and intensity will keep on till you stop the cycle function by pressing the cycle icon again. The default cycle countdown time is 20 minutes, the cycle function will make the countdown time repeat again and again. You can also reset the countdown time by yourself.
- If you enjoy certain mode and don't want to make any change further, you can press the ON/OFF button once slightly, the screen will be locked and the backlight is off in order to save more energy, Or press the lock icon on the screen, the screen will be locked but the backlight is still on. In locking status, the unit is still on with intensity output but all function icons are locked, you are not

able to change the mode, intensity etc. It will avoid screen touching by mistake. To unlock the screen, you can press the ON/OFF once slightly or press the lock icon on the screen.

- Pressing the countdown time on the screen, you can set the timer from 10 minutes to 60 minutes.
- When power bar in the battery icon disappear, the unit is going to run out power, it should be charged as soon.
- When you turn on the unit, if the unit test no load, for example you didn't connect the unit well with body, you are not able to adjust the intensity.

**Recommended practice:**

- Duration suggested for each skin area is 20 min and 2 times per day. Consult with your physician for longer and more frequent uses.
- Start from the lowest intensity and gradually adjust the intensity to a comfortable level at a scale from 1 to 20.
- Good skin care is important for a comfortable use of device. Be sure the treatment site is clean of dirt and body lotion.
- Keeping the electrode in the storage bag after use will extend its lifespan. The electrode is disposable and should be replaced when it loses the adhesiveness. To purchase additional electrodes, please contact the seller.

## **Cleaning and maintenance**

---

Please use wipe to clean the device first, and then use the dry cloth to wipe it again. The electrodes coming with the device are disposable, and should be replaced when their adhesiveness becomes worse. Contact the seller for replacements.

## Trouble Shooting

If your device is not operating properly, please check below for common problems and suggested solutions. If the recommended action does not solve the problem, please contact the seller.

Problem	Possible Cause	Solution
The intensity is not felt with a very weak intensity level	Pads are not attached to the body firmly	Attach both pads firmly to the skin
	The transparent films are still stuck to the pads	Peel off film on the adhesive surface of pads
	The pads stack together or overlap	Do not stack pads together or overlap pads
	The intensity setting is getting weak	Increase the intensity level
	The battery capacity is low	Change the battery
The skin turns red or the skin feels irritated	The adhesive surface of the pads is dirty or dry	Wash adhesive surface of pads gently with your fingertips for about 3 seconds under slow running water
	The therapy time is too long or the intensity is set too high	Reduce the application time or reduce the intensity
	The electrode pad surface if worn out	Replace electrode pad
No power source	The battery capacity is depleted	Change the battery
Power cuts off during use	The battery is weak	Change the battery
It is difficult to attach the pad to the skin	Have you remove the transparent film from the pad?	Peel off film on the adhesive surface of pads
	Was the pad applied immediately after washing?	Dry the pad
	Is the adhesive surface of the pad damaged?	Replace the pad
Adhesive surface of pad is not sticky	Pads get deteriorative	Contact the vendor for replacements.
	Were the pads stored under high temperature, high humidity, or sunshine	Replace the pad.

## Recommended Use Positions

---

SHOULDERS

BACK

BUTTOCKS

ABDOMEN


ARMS

LEGS

FEET

## Contact Information

---

**JKH Health Co., Ltd.**

Address: 4-5F, Building 12, Hengmingzhu Industrial Park, Xinqiao Tongfuyu Industrial Area, Shajing, Baoan District, Shenzhen, China

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